

FIRST PROFESSIONAL SPOT TEST

B.U.K LEVEL 400 2016/17

NOTE


The time for each spot identification is 30seconds

SPOT A = CUVETTE

Path length is defined as the distance that light travels through a sample in an analytical cell.



Cuvette is used to hold samples for spectroscopic analysis




SPOT B = WAVELENGTH SELECTOR FILTER



Used in the colorimeter to select wavelength which the solute absorbs the most, in order to maximize accuracy

SPOT C = WESTEGREN STAND WITH WESTEGREN TUBE

It is used for ESR estimation



SPOT D

BLUE TOP = lithium Heparin container

YELLOW TOP = Fluoride oxalate container



- **Sodium heparin** = For serum/chemistry (urea and electrolytes)
- **Fluoride oxalate** = Use for glucose estimation

SPOT E = HEMOGLOBIN PIPETTE WITH TUBING



SPOT F = WIRE LOOP



used to inoculate test samples into culture media for bacterial or fungal culture (for streaking inoculum unto agar plates)

SPOT G = LABORATORY INCUBATOR



are mainly used for growing microbial cultures

SPOT H = ASCARIS LUMBRICOIDES



Definitive host is man

SPOT I = BOUINS FLUID

- It is **composed** of picric acid, acetic acid and formaldehyde in an aqueous solution.
- It is especially good for gastrointestinal tract biopsies because this fixative allows crisper and better nuclear staining than 10% neutral-buffered formalin
- It is prepared as follows:
 - Picric acid, saturated aqueous solution - 75ml;
 - Formalin, 40% aqueous solution - 25ml;
 - Glacial Acetic acid, 1 - 5ml.

SPOT J = SAPA

(Saturated Alcoholic Picric Acid)

Constituents

- Picric acid
- Alcohol (e.g Ethanol)

NOTE

- Picric acid solution is a saturated alcoholic formulation for use as a background stain in where it will produce a yellow background or as a fixative, it is often used in conjunction with **Trichrome stains** because its use as fixative or in the staining protocol can very much enhance the stain picture obtained from the trichrome stain (stains erythrocyte).

OTHER POSSIBLE SPOTS

SPECTROPHOTOMETER



It is used to measure the intensity of light absorbed after it passes through a sample solution

CUVETTE HOLDER



It is designed to accommodate(hold) standard micro and macro cuvettes within the spectrophotometer.

COLORIMETER



It is used for measuring transmittance and absorbance of light passing through a sample solution

FLAME PHOTOMETER



It is used to determine the concentration of certain metal ions – Na, K, Cl & Ca

HOT AIR OVEN




Hot air oven is usually used for heating, drying & sterilizing items

AUTOCLAVE



It is used to sterilize equipment and supplies by subjecting them to high-pressure saturated steam at 121°C for 15-20 min depending on the size of the load and the contents.

HEMATOCRIT CENTRIFUGE



It is used for determination of volume fractions of erythrocytes in blood (PCV) and also for separation of micro blood and solutions

LABORATORY CENTRIFUGE



It is used to separate blood components

DIGITAL WATER BATH




- It is use to incubate samples in water at constant temperature over a long period of time

ELECTROPHORESIS MACHINE



- It is used to separate macromolecules either nucleic acids or proteins on the basis of size, electric charge and other physical properties.

ELECTROPHORESIS TANK



- It is where the electrophoresis buffer is poured.

ELECTROPHORESIS MACHINE WITH ELECTROPHORESIS TANK ON TOP



WEIGHING MACHINE



- It is used to measure the weight of an item.

HOT PLATE



- It is for Heat drying glass slides (or other glass wares)

MAGNETIC STIRRER WITH BEAKER ON TOP




It is used for stirring heated liquid automatically

HOT PLATE AND STIRRER



It is used for heating and automatic stirring

MICROSCOPE



OBJECTIVE LENS



- It relays a real image of the object to the eyepiece

EYE PIECE



For magnifying the primary image produced by the objective lens

4X OBJECTIVE LENS



For scanning parasites

10x OBJECTIVE LENS (Yellow)



For focusing wet preparations

40x OBJECTIVES (Blue)



For viewing wet slides and focusing stained slides

100x OBJECTIVE LENS (white)



For viewing stained slides

CONDENSER



For supplying the object with a cone of light correct size and formular

IRIS DIAPHRAGM



For regulating the angle of cone of light entering the condenser

PLAIN CONTAINER



For storing blood samples for serology

LITHIUM HEPARIN CONTAINER



For storing blood for renal function test

FLUORIDE OXALATE CONTAINER



For storing blood for glucose test

EDTA CONTAINER



For storing blood samples for lab test

VACUTAINERS

PLASTIC TUBE TOP COLOR AND ORDER OF DRAW – I

Draw Blood Culture bottles first, then proceed with blood tube order of draw (1-12) below:			
1. 	LIGHT BLUE top plastic tube PT, PTT, Fibrinogen, Fibrin D-Dimer, other Coagulation Testing Note: invert gently 3 - 4 times	7. 	BRIGHT GREEN top (SODIUM HEPARIN) plastic non-gel tube Mycobacteriology (AFB) Blood Culture, HLA-B27, Chromosome Studies
2. 	GOLD gel plastic tube Most Chemistry tests & Immunology Tests, Hepatitis Tests, Serologies (Do not use for Troponin, BNP)	8. 	LAVENDER top plastic tube Hematology: CBC, Platelet, Sed. Rate Chemistry: CD4, CD8, G6PD, Hemoglobin A1C & Hemoglobin Variants
3. 	RED top plastic tube For tests requiring serum Note: contains clot activator	9. 	WHITE top plastic tube (PPT) Hepatitis and HIV Viral Loads, BNP
4. 	ROYAL BLUE top plastic tube Copper, Zinc, Trace Elements	10. 	PINK top plastic tube for Blood Bank <u>ONLY</u>.
5. 	LIGHT GREEN top (LITHIUM HEPARIN) gel plastic tube Troponin, Metabolic Panels, Lipid, Liver Panels, Ammonia (ice), HIV Rapid Anti-body	11. 	TAN top plastic tube Lead
6. 	DARK GREEN top (LITHIUM HEPARIN) plastic non-gel tube Ionized Calcium (not part of blood gas), Ammonia (ice)	12. 	GRAY top plastic tube Glucose, Lactate (Lactic Acid) on ice
* IMPORTANT *: Please follow the correct order of draw as numbered above and thoroughly mix all specimens (except Light Blue top) by inversion 8 - 10 times.			

PLAIN VACUITAINER




It is used in serum/chemistry and serology tests

GEL ACTIVATOR



It is used in serum/chemistry, immunology, Hepatitis test

THROMBIN




It is used for stat serum/chemistry testing

EDTA VACUITAINER (PURPLE)



It is used in plasma/ Haematology -FBC

EDTA VACUITAINER (PINK)



SODIUM/LITHIUM HEPARIN VACUITANER



It is used in plasma/chemistry tests –urea & electrolytes

ROYAL BLUE TOP VACUITAINER



It is used for trace metal elements(-Cu, Zn) analysis & Toxicology

TAN



It is used for lead determinations.

SODIUM CITRATE VACUITANER



It is used in plasma/coagulation tests

ESR VACUITAINER



It is used in plasma/ ESR –Haematology tests

FLUORIDE OXALATE VACUITANER



It is used in plasma/glucose testing/Lactic acid, blood alcohol

LIGHT YELLOW VACUITAINER




Acid-citrate- dextrose (ACD): It is Used for used blood bank studies, HLA phenotyping, and paternity testing . Or SPS in serum microbiology/culture

WHITE TOP VACUITAINER




It is used for hepatitis and HIV viral load, BNP

SLIDE STAINING RACK




- It is used for the simultaneous **staining** of microscope slides

STAINING RACK



- It is used for the simultaneous **staining** microscope slides

DRYING RACK (WITH SLIDES)



- It is used for draining slides in a laboratory

STAINING DISH



staining dishes are **used** to **stain** cells and tissues on microscope slides.

GLASS STAINING JAR WITH DIP STAINING RACK (SIDE)



DIP STAINING RACK



- It is used for the simultaneous **staining** microscope slides by dipping vertically in a staining dish.

COPLIN JAR



- **Staining jars** are used to **stain** cells and tissues on microscope slides

TEST TUBE



- *Test tubes* are used to hold, mix, and heat chemical experiments. They are used as homes for microorganisms when people want to culture (grow) them.

CENTRIFUGE TUBE




For spinning urine samples

WESTEGREN TUBE ON A CUP



For ESR estimation

WESTEGREN CUP AND TUBE



For ESR estimation

DURHAM TUBE



For detecting gas production in sugar fermentation media by microorganisms
(the tube is placed in an inverted fashion so that gases produced get trapped in it and do not float away to the surface)

BIJOU BOTTLE



For carrying out biochemical tests

McCartney BOTTLE



McCartney is used for simultaneous solid and liquid cultures.

TEST TUBE BRUSH



For washing test tubes

TEST TUBE HOLDER




For carrying test tubes

TEST TUBE HOLDER WITH TEST TUBE




TEST TUBE RACK (plastic)



It is used to **hold** upright **multiple** test tubes at the same time

TEST TUBE RACK (WOODEN)



It is used to **hold** upright **multiple** test tubes at the same time

TEST TUBE RACK (Metal)



It is used to **hold** upright **multiple** test tubes at the same time

AUTOMATIC PIPETTE



used in transferring measured volume of liquid or solution

MULTI CHANNEL PIPETTE



It is used in pipetting various samples at a time which increases output

AUTOMATIC PIPETTE TIP



It is used to avoids cross contamination of pipettes

PASTEUR PIPETTE



used in transferring small volume of liquid or solution

RBC PIPETTE



use for dilution of blood with RBC diluting fluid to count the number of red blood cells

WBC PIPETTE



use for dilution of blood with WBC diluting fluid to count the number of white blood cells

DIFFERENCES BETWEEN RBC AND WBC PIPETTE

	RBC pipette	WBC pipette
1)	It has a red bead	It has a white bead
2)	It has graduations upto mark 101	It has graduations upto mark 11
3)	Size of bulb is larger	Size of bulb is smaller
4)	Size of lumen is smaller	Size of lumen is larger

GLASS PIPETTE



For pipetting and dispensing a variable volume of aqueous liquid. E.g during titration

PIPETTE FILLER



it is Use to safely fill glass and plastic pipettes.

PIPETTE FILLER



- it is Use to safely fill glass and plastic pipettes.


TRIPOD STAND



For placing objects for heating


WIRE GAUZE

Wire gauze with filter



- It is used to support the beakers or other glassware or flasks during heating

• Wire gauze without filter




TRIPOD STAND WITH WIRE GAUZE



- It is used to support or hold the flasks and beakers during heating and the wire gauze protects them from direct heating

ANAEROBIC JAR



For culturing of anaerobic organism




Figure 6.6

WASH BOTTLE



For rinsing test tubes, cuvettes

TISSUE FLOATING BATH



It is used to stretch sections and removes wrinkles & folds

VORTEX MIXER



It is used to mix small vials of liquid

BLOOD MIXER



It is used for complete mixing of blood prior to testing

NEUBAEUR CHAMBER



For manual counting of Wbc, Rbc, platelets and semen analysis on the microscope

COMPARATOR



- It is used for matching the colour of acid haematin formed in the haemoglobinometer tube.

SAHLI'S HEMOGLOBINOMETER

- Or Haemometer



L-R:
Pipette
Sahli's Standard
Hemometer Tube
Stirring Rod
Dropper

- It is **used** to determine the hemoglobin content of the blood

ANALYTICAL BALANCE



It is used for accurate weighing of samples and precipitates mass in the sub-milligram range

DISTILLER



It is used to produce high quality water through distillation

DEIONIZER



It is used to remove mineral salts from water by ion exchange resins

MICROTOME KNIFE



It is used to cut thin sections of tissue and trimming



CANNISTER



It is used for culturing anaerobic organisms/sterilizing and drying plates in autoclave or hot air oven